DUROLANE INSTRUCTIONS FOR USE

Read manufacturer **INSTRUCTIONS FOR USE** completely and follow the instructions prior to use.

Contents		
Each ml contains:		
Hyaluronic acid stabilized	20 mg	
Phys sodium chloride solution, pH7	q.s.	

Indications

Treatment of patients with knee osteoarthritis pain, restricted to conservative non-drug treatment and general analgesics such as Acetaminophen use invalid.

Description

DUROLANE contains 20 mg/ml of stabilized non-animal hyaluronic acid in buffered physiological sodium chloride solution pH7. DUROLANE is a sterile, transparent viscoelastic gel supplied in a 3 ml glass syringe. The product is for single use only.

Hyaluronic acid belongs to a group of very few substances which are identical in all living organisms. It is a natural poly-saccharide that is present throughout the tissues of the body, with particularly high concentrations in the synovial fluid and the skin. DUROLANE is composed of biosynthetically produced hyaluronic acid which has been purified and stabilized. DUROLANE is degraded in the body by the same metabolic pathway as endogenous hyaluronic acid.

Mode of Action
The body's hyaluronic acid constitutes a natural part of the synovial fluid and acts in the joints both as a lubricant of cartilage and ligaments and as a shock absorber. It is known that the synovial fluid in joints affected by osteoarthritis has a much lower viscosity and elasticity than in healthy joints. Injections of hyaluronic acid in the joint to restore the viscosity and elasticity can diminish the pain and improve the mobility of the joint.

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Dosage
DUROLANE is a single dose preparation and should only be injected once per treatment course. The recommended dose is 3 ml (one syringe) per knee joint.

Contraindications
None known.

Warnings

- DUROLANE should not be injected if the knee joint is infected or severely inflamed.
- DUROLANE should not be injected if there is an active skin disease or infection present at or near the injection site.
- DUROLANE should not be injected intravascularly or extra-articularly or in the synovial tissues or capsule.

Precautions

- DUROLANE should be used with caution in patients with venous or lymphatic stasis present in the leg.
- DUROLANE has not been tested in pregnant or lactating women or in children.
- If treatment is bilateral, a separate syringe of DUROLANE must be used for each knee.

- As with any invasive joint procedure there is a small risk of infection when injecting DUROLANE.

- DUROLANE should not be injected if the patient is known to be sensitive to hyaluronic acid based products.

- Local anaesthetics should not be used if the patient is known to be allergic or sensitive to local anaesthetics.

- In clinical studies, reinjections in the knee have not been studied with a shorter interval between first and second injection than 6 months.

- Clinical use beyond two doses of DUROLANE has not been studied.

- DUROLANE should be used with caution in patients with pre-existing chondrocalcinosis as injection may lead to an acute attack of the condition.

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Adverse Events
The majority of the reported adverse reactions in clinical studies were described as transient pain, swelling and/or stiffness–localized to the knee and were of mild or moderate intensity with a median duration of one week.

In a few patients symptoms of pain and/or swelling/stiffness localized to the knee lasted for more than 3 weeks but in these cases the observed symptoms were not distinguishable from fluctuations in the underlying osteoarthritis condition.

Adverse events must be reported to the local Bioventus representative.

Interactions
The safety and effectiveness of DUROLANE concomitantly with other intra-articular injectables have not been established.

Administration

General administration information

- DUROLANE should only be injected by an authorized physician (or in accordance with local legislation) in facilities well suited for intra-articular injections.

- DUROLANE should be injected using strict aseptic technique.

- The injection site should be swabbed with alcohol or other suitable antiseptic solution before injection.

- The route for injection with DUROLANE should be chosen so that damage to adjacent vital structures is avoided.

- DUROLANE should be injected into the joint cavity only.

- Remove joint effusion, if present, before injecting DUROLANE. The same needle should be used for both removal of effusion and injection of DUROLANE.

- The recommended needle size is 18 to 22 G and with adequate length.

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Please inform your patient that:

- As with any invasive joint procedure it is recommended to avoid strenuous activity (e.g. tennis, jogging or long walks) the first two days after the injection.

- Some transient reactions related to the injection of DUROLANE, such as pain and/or swelling/stiffness of mild to moderate intensity during the first week following the injection can be anticipated. If the symptoms last for more than a week a physician should be contacted.

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Performance
• Clinical data indicate significant mean benefits, such as improvement in knee pain and physical function, versus baseline at 6 months post-treatment. Patients were retreated after 6 months. Retreatment did not give rise to an increased rate of adverse events.

- The half life of the product in human knees is approximately four (4) weeks, meaning that about one-fourth of the dose is left after 8 weeks.

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How Supplied
DUROLANE is supplied in a 3 ml glass syringe with a Luer-lok fitting, packed in a blister pack. The contents of the syringe, stabilized non-animal hyaluronic acid gel, are sterile.

DUROLANE is intended for single use and should not be re-sterilized. It should be used immediately after the syringe has been removed from its packaging. If the blister package or syringe is opened or damaged, do not use.

Shelf life and Storage
DUROLANE should be stored, in its original packaging, up to 30 °C. The expiry date is indicated on the package. Protect from freezing.

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Manufacturing site
Q-MedAB, Seminariegatan, 21, SE-752 28 Uppsala, Sweden

For
Bioventus LLC, 4721 Emperor Blvd., Suite 100, Durham, NC 27703 USA

Company name: Emergo Taiwan Ltd. located at 10F., No.225, Sec. 3, Beixin Rd., Xindian Dist., New Taipei City 231, Taiwan

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IF THE PACKAGE IS DAMAGED, DO NOT USE

*DUROLANE and NASHA are registered trademarks of Galderma S.A.

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Bioventus and the Bioventus logo are trademarks of Bioventus LLC.

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